




Pedsan

GENERAL CATALOGUE NO. 6





For more than forty years on the market, we have always had as our primary target the customer satisfaction. Our work, in fact, is focused on a careful selection of raw materials, which represent the highest quality standards and a sustainable and extremely flexible production system, able to better meet the changing needs.




Our strength is the artisan production, which can count on the experience of skilled operators and that is entirely carried out in our factory in Perugia, as the real Made in Italy asks. Our commitment and your confidence will allow us to successfully face new challenges and to develop interesting collaboration with you.


■ FLAT / PRONATED FOOT INSOLES	P. 2 - 17
■ FOREFOOT PATHOLOGIES INSOLES	P. 18 - 34
■ PLANTAR FASCIITIS AND CALCANEAL SPUR INSOLE	P. 35 - 38
■ RHEUMATIC AND DIABETIC FOOT INSOLES	P. 39 - 44
■ SPORTS	P. 45 - 50
■ VARIOUS PATHOLOGIES INSOLES	P. 51 - 57
■ PHYSIOLOGICAL INSOLES	P. 58 - 63
■ MONOBLOCK INSOLES	P. 64 - 65
■ ORTHOKIT / SILIPED	P. 66
■ DYSMETRIA WEDGE	P. 67
■ CORRECTIVE ELEMENTS	P. 68 - 69
■ FOOTPRINT / GLUING SYSTEM	P. 70
■ SHAPED COVERING / INSOLE BAGS / MINIROLL	P. 71
■ BLOCK FOR CAD/CAM INSOLES	P. 72 - 74
■ RAW MATERIALS	P. 75 - 87

ART. B550MHEVA

Recommended for any flat-footedness degrees treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base


SIZES: from 20 to 47

ART. B551HEVA

Recommended for any flat-footedness degrees treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base


SIZES: from 20 to 47

ART. A551HEVA

Recommended for any flat-footedness degrees with mild metatarsalgia treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base

SIZES: from 20 to 47

ART. A550MHE

Recommended for any flat-footedness degrees with mild metatarsalgia treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base

SIZES: from 20 to 47

ART. A550VHE

Recommended for any flat-footedness degrees with mild metatarsalgia treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base

SIZES: from 20 to 46

ART. B300MHE MST VARIANT

Recommended for valgus flat foot treatment



AVAILABLE COLOURS:



Notes: _____


SIZES: from 20 to 46




SIZES: from 20 to 46

ART. B300MHEVA

Recommended for valgus flat foot treatment



AVAILABLE COLOURS:




Notes: _____




SIZES: from 22 to 48

ART. B553HEVA

Recommended for the pronated flat foot treatment




AVAILABLE COLOURS:



Notes: _____

ART. A500MHE

Recommended for any flat-footedness degrees with mild metatarsalgia treatment



AVAILABLE COLOURS:




Notes: _____


SIZES: from 20 to 46

ART. B500MHEVA

Recommended for any flat-footedness degrees treatment



AVAILABLE COLOURS:



Notes: _____


SIZES: from 20 to 46



SIZES: from 20 to 46

ART. B500VHE

Recommended for any flat-footedness degrees treatment



COLORI DISPONIBILI:




Notes: _____




SIZES: from 20 to 46

ART. B502HEVA

Recommended for any flat-footedness degrees treatment




AVAILABLE COLOURS:



Notes: _____

ART. B600HE

Recommended for pronated flat foot treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base

SIZES: from 20 to 46

ART. 90HE

Recommended for pronated flat foot treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base

SIZES: from 20 to 46

ART. A90HE

Recommended for pronated flat foot with mild metatarsalgia treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base


SIZES: from 20 to 46

ART. 90HEVA

Recommended for pronated flat foot treatment



AVAILABLE COLOURS:



Notes: also available with prefinished base


SIZES: from 20 to 46



SIZES: from 20 to 46

ART. 95HEVA

Recommended for pronated flat foot treatment




Notes: _____



SIZES: from 20 to 46

ART. 9055H

Recommended for pronated flat foot treatment




Notes: _____

ART. B700E MCC VARIANT

Recommended for pronated flat foot treatment

Microcel top layer
Thickness: 5 mm
Density: approx. 300 Kg/m³
Hardness: approx. 50 shore A

Microcel three-quarter base
Thickness: 10 mm
Density: approx. 300 Kg/m³
Hardness: approx. 50 shore A



AVAILABLE COLOURS:



Notes: _____


SIZES: from 22 to 48

ART. B700E GRADRC VARIANT

Recommended for pronated flat foot treatment

EVA Granite top layer
Thickness: 5 mm
Density: approx. 160 Kg/m³
Hardness: approx. 30 shore A

Durocel three-quarter base
Thickness: 10 mm
Density: approx. 320 Kg/m³
Hardness: approx. 60 shore A



AVAILABLE COLOURS:



Notes: _____

SIZES: from 22 to 48

ART. GP

Recommended for serious pronated flat foot treatment

Polyethylene layer
Density: 950 kg/m³
Hardness: 62 shore D



Heel cup with
wraparound edges

AVAILABLE COLOURS:



Notes: also available semifinished (with lateral edge not worked)

SIZES: from 24 to 46 (only even numbers)

ART. QSHE

Recommended for navicular drop treatment

Micro-perforated Veolene
top layer
Thickness: 4 mm
Density: 140 Kg/m³
Hardness: 29 Shore A

Latex foam quarter of sphere
Density: 300 Kg/m³
Hardness: 25 Shore E

Duroform 1,4 mm
reinforcement from
the heel to metatarsals



External
reinforcement

Microcel three-quarter base
Thickness: 10 mm
Density: approx. 300 Kg/m³
Hardness: approx. 50 Shore A



AVAILABLE COLOURS:



Notes: also available with prefinished base


SIZES: from 20 to 40



SIZES: from 20 to 43

ART. LHE

Recommended for knees and rearfoot valgus treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base



SIZES: from 23 to 45

ART. LCG

Recommended for pronated flat foot treatment



AVAILABLE COLOURS:



Notes:

ART. TMA550

Recommended for any flat-footedness degrees with mild metatarsalgia treatment

Cork Mousse top layer
Thickness: 3 mm
Density: approx. 420 Kg/m³
Hardness: approx. 28 Shore A

Thermoformed drop shaped metatarsal pad

Thermo mousse three-quarter base
Thickness: 10mm
Density: 490 kg/m³

AVAILABLE COLOURS:



Notes: also available with prefinished base

SIZES: from 20 to 47

ART. TMB550C

Recommended for any flat-footedness degrees treatment

Vegetable-tanned leather entire base
First choice leather
Thickness: 1,8 mm

Cork Mousse top layer
Thickness: 3 mm
Density: approx. 420 Kg/m³
Hardness: approx. 28 Shore A

Thermo mousse inner layer
Thickness: 10 mm
Density: 490 Kg/m³

AVAILABLE COLOURS:




Notes: _____

SIZES: from 20 to 47

ART. A2M

Recommended for any flat-footedness degrees with mild metatarsalgia treatment



SIZES: from 20 to 47


AVAILABLE COLOURS:



Notes: _____

ART. B2M

Recommended for any flat-footedness degrees treatment



SIZES: from 20 to 47


AVAILABLE COLOURS:



Notes: _____

ART. A3M

Recommended for any flat-footedness degrees with mild metatarsalgia treatment



SIZES: from 20 to 46


AVAILABLE COLOURS:



Notes: _____

ART. B3M

Recommended for any flat-footedness degrees treatment



SIZES: from 20 to 46

AVAILABLE COLOURS:




Notes: _____



SIZES: from 20 to 43

ART. L

Recommended for knees and rearfoot valgus treatment



AVAILABLE COLOURS:




Notes: _____



SIZES: from 20 to 46

ART. PA

Recommended for pronated flat foot treatment




AVAILABLE COLOURS:



Notes: _____

ART. P48005

Recommended for claw foot with metatarsalgia treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base


SIZES: from 34 to 48

ART. P48006

Recommended for claw foot with metatarsalgia treatment



AVAILABLE COLOURS:




Notes: _____

SIZES: from 34 to 48

ART. PED06B-R

Recommended for claw foot with metatarsalgia treatment



AVAILABLE COLOURS:




Notes: _____

SIZES: from 34 to 48

ART. PED08B-R

Recommended for claw foot with metatarsalgia and calcaneodynia treatment



AVAILABLE COLOURS:



Notes: _____


SIZES: from 34 to 48



SIZES: from 34 to 48

ART. PED02UB

Recommended for claw foot with metatarsalgia and calcaneodynia treatment



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 48

ART. D779THE

Recommended for claw foot with metatarsalgia treatment




AVAILABLE COLOURS:



Notes: _____

ART. 33FCEA

Recommended for claw foot with metatarsalgia treatment



AVAILABLE COLOURS:




Notes: _____

SIZES: from 34 to 48

ART. PED05B

Recommended for claw foot with metatarsalgia and plantar fasciitis treatment



AVAILABLE COLOURS:



Notes: _____


SIZES: from 34 to 48



SIZES: from 34 to 42

ART. D771HE

Recommended for claw foot with metatarsalgia treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base



SIZES: from 34 to 48

ART. U881HGRA

Recommended for claw foot with metatarsalgia treatment




AVAILABLE COLOURS:



Notes: also available with prefinished base

ART. U881EBV

Recommended for claw foot with metatarsalgia treatment



AVAILABLE COLOURS:




Notes: _____


SIZES: from 34 to 48

ART. U881HEVA

Recommended for claw foot with metatarsalgia treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base

SIZES: from 34 to 48

ART. U882HE

Recommended for claw foot with metatarsalgia treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base

SIZES: from 34 to 47

ART. U800MED

Recommended for claw foot with metatarsalgia treatment



AVAILABLE COLOURS:



Notes: also available with prefinished base


SIZES: from 34 to 46



SIZES: from 34 to 47

ART. P33004

Recommended for claw foot with metatarsalgia treatment



AVAILABLE COLOURS:



LIME


Notes: also available with prefinished base



SIZES: from 34 to 47

ART. P33005

Recommended for claw foot with metatarsalgia treatment



AVAILABLE COLOURS:



BM


Notes: also available with prefinished base



SIZES: from 34 to 49

ART. NR01

Recommended for claw foot with Morton's Neuroma treatment



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 48

ART. NR02

Recommended for claw foot with metatarsalgia and Morton's Neuroma treatment



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 42

ART. D749HX

Recommended for valgus claw foot with metatarsalgia treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base



SIZES: from 34 to 47

ART. NATURE

Recommended for claw foot with metatarsalgia treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base



SIZES: from 34 to 42

ART. D748HE

Recommended for claw foot with metatarsalgia and hammer toes treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base



SIZES: from 34 to 48

ART. U848HE

Recommended for claw foot with metatarsalgia and hammer toes treatment




AVAILABLE COLOURS:



Notes: also available with prefinished base

ART. U881HK

Recommended for claw foot with metatarsalgia treatment



SIZES: from 34 to 48


AVAILABLE COLOURS:



Notes: also available with prefinished base

ART. U881HEA

Recommended for claw foot with metatarsalgia treatment



SIZES: from 34 to 48

AVAILABLE COLOURS:



Notes: also available with prefinished base

ART. D700HEVA

Recommended for claw foot with metatarsalgia of second, third and fourth metatarsal treatment


EVA top layer
Thickness: 4 mm
Density: approx. 200 Kg/m³
Hardness: 35/45 Shore A

Drop metatarsal pad with a latex foam insert

Sintex 1,4 mm reinforcement from heel to metatarsals

Stiff Veolene three-quarter base
Thickness: 10 mm
Density: 170 Kg/m³
Hardness: 40 Shore A

AVAILABLE COLOURS:



Notes: also available with prefinished base

SIZES: from 34 to 42

ART. U800HE

Recommended for claw foot with metatarsalgia of second, third and fourth metatarsal treatment

Micro-perforated Veolene top layer
Thickness: 4 mm
Density: 140 Kg/m³
Hardness: 29 Shore A

Drop metatarsal pad with a latex foam insert

Sintex 1,4 mm reinforcement from heel to metatarsals

Stiff Veolene three-quarter base
Thickness: 10 mm
Density: 170 Kg/m³
Hardness: 40 Shore a

AVAILABLE COLOURS:




Notes: also available with prefinished base

SIZES: from 34 to 46

ART. D110E

Recommended for claw foot with metatarsalgia of second, third and fourth metatarsal treatment



SIZES: from 34 to 42


AVAILABLE COLOURS:



Notes: also available with prefinished base

ART. D117E

Recommended for claw foot with metatarsalgia treatment



SIZES: from 34 to 42

AVAILABLE COLOURS:




Notes: also available with prefinished base



SIZES: from 34 to 42

ART. D1/LU

Recommended for claw foot with metatarsalgia of second, third and fourth metatarsal treatment



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 46

ART. U1

Recommended for claw foot with metatarsalgia of second, third and fourth metatarsal treatment



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 42

ART. D71/N

Recommended for metatarsalgia treatment



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 42

ART. D76/N

Recommended for metatarsalgia of second, third and fourth metatarsal and Morton's Neuroma treatment



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 48

ART. P33006

Recommended for claw foot with metatarsalgia treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base




SIZES: from 34 to 48

ART. F848EVA

Recommended for calcaneal spur, plantar fasciitis, metatarsalgia and hammer toes treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base



SIZES: from 34 to 48

ART. F848E

Recommended for calcaneal spur, plantar fasciitis, metatarsalgia and hammer toes treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base



SIZES: from 34 to 48

ART. F848E-PLUS

Recommended for calcaneal spur, plantar fasciitis and metatarsalgia treatment



AVAILABLE COLOURS:




Notes: _____




SIZES: from 34 to 49

ART. P52003

Recommended for claw foot with calcaneal spur treatment



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 50

ART. F852E

Recommended for plantar fasciitis and calcaneal spur treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base




SIZES: from 34 to 48

ART. PED05UB-R

Recommended for calw foot with plantar fasciitis and metatarsalgia treatment



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 48

ART. P33003

Recommended for claw foot with metatarsalgia and calcaneal spur treatment



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 46

ART. F100E

Recommended for calcaneal spur treatment



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 48

ART. RD01

Recommended for rheumatic and diabetic foot pathologies prevention and treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base



SIZES: from 34 to 50

ART. RD02

Recommended for rheumatic and diabetic foot pathologies prevention and treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base



SIZES: from 34 to 48

ART. RD03

Recommended for rheumatic and diabetic foot pathologies prevention and treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base



SIZES: from 34 to 48

ART. RD04

Recommended for rheumatic and diabetic foot pathologies prevention and treatment



AVAILABLE COLOURS:




Notes:



SIZES: from 34 to 46

ART. P46001 TOT VARIANT

Recommended for rheumatic and diabetic foot pathologies prevention and treatment



AVAILABLE COLOURS:




Notes: also available with a three-quarter base



SIZES: from 34 to 46

ART. P46002

Recommended for rheumatic and diabetic foot pathologies prevention and treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base



SIZES: from 34 to 50

ART. 900EMTE

Recommended for rheumatic and diabetic foot pathologies prevention and treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base



SIZES: from 34 to 48

ART. D771TMTE

Recommended for rheumatic and diabetic foot pathologies prevention and treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base



SIZES: from 34 to 48

ART. RD05

Recommended for rheumatic and diabetic foot pathologies prevention and treatment



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 48

ART. DINATEK UB HVNCX VARIANT

Recommended for rheumatic and diabetic foot pathologies prevention and treatment



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 50

ART. DINATEK U

Recommended for rheumatic and diabetic foot pathologies prevention and treatment



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 48

ART. DINATEK TB MCR VARIANT

Recommended for rheumatic and diabetic foot pathologies prevention and treatment



AVAILABLE COLOURS:




Notes: _____




SIZES: from 34 to 46

ART. PLSPORT 01

Recommended for backfoot instability with metatarsalgia treatment



AVAILABLE COLOURS:




Notes: _____




SIZES: from 34 to 46

ART. PLSPORT 02

Recommended for claw foot with metatarsalgia treatment



AVAILABLE COLOURS:




Notes: _____




SIZES: from 34 to 48

ART. PLSPORT 03

Recommended for claw foot with metatarsalgia and talalgia treatment



AVAILABLE COLOURS:




Notes: _____




SIZES: from 34 to 48

ART. PLSPORT 04

Recommended for claw foot with metatarsalgia and talalgia treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base _____



SIZES: from 34 to 48

ART. PLSPORT 06

Recommended for claw foot treatment with an amortised insert for calcaneal, cuboid bone and metatarsals



AVAILABLE COLOURS:




Notes: also available with prefinished base



SIZES: from 34 to 50

ART. PLSPORT 07

Recommended for plantar fasciitis and calcaneal spur treatment



AVAILABLE COLOURS:




Notes: also available with prefinished base



SIZES: from 34 to 46

ART. RUNNING 01

Recommended for running, let a total load distribution with a shock-absorbing insert till metatarsals



AVAILABLE COLOURS:




Notes: _____




SIZES: from 34 to 46

ART. RUNNING 02

Recommended for heel strike and Achilles tendinopathy of the runner, let a total load distribution with a full shock-absorber base



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 46

ART. 201NE

Recommended for valgus claw foot with mild metatarsalgia treatment, with a shock absorber insert till metatarsals



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 47

ART. ORTHOLIFE

Recommended for backfoot and forefoot pathologies of athlete prevention



AVAILABLE COLOURS:



Notes: _____



SIZES: from 34 to 47

ART. ORTHOSPORT

Recommended for backfoot and forefoot pathologies of athlete prevention with a total shock absorber base in Noene®

Evaform 180 top layer
CAD / CAM worked
Density: approx. 180 Kg/m³
Hardness: approx. 34 Shore A

Metatarsal bar

Noene® total shock absorber base
Thickness: 2 mm
Density: approx. 640 Kg/m³
Hardness: approx. 34 Shore A

Heel cup

AVAILABLE COLOURS:



Notes: also available in LATEX FOAM, Black EVAFORM 100, Grey EVAFORM 120, Beige EVAFORM 160, Brown EVAFORM 240, Fleshtone VEOLENE, MULTISTRATOS, Marbled EVACOLOR and MICROCEL



SIZES: from 35 to 48

ART. PEDSPORT

Recommended for a backfoot stabilization with a propulsive action

Airlastic, visco-elastic material,
top layer with a microfibre cover
Thickness: 3 mm

Thermoformed drop-shaped
metatarsal pad

Millenium and Duroform
calcaneal stabilizer

Duroform double
reinforcement
layer

Suede Microfibre lower cover
Thickness: 0,8 mm

Poron® XRD shock
absorber insert
Thickness: 3 mm
Density: 240 Kg/m³

AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 48

ART. PLCAD CE VARIANT

High thickness CAD / CAM worked insole. Adjustable to various pathologies



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 48

ART. PLCAD MST VARIANT

High thickness CAD / CAM worked insole. Adjustable to various pathologies



AVAILABLE COLOURS:



Notes: _____



SIZES: from 34 to 48

ART. PLCAD EVAMIX VARIANT

High thickness CAD / CAM worked insole. Adjustable to various pathologies.



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 48

ART. PLCAD CEGR VARIANT

High thickness CAD / CAM worked insole. Adjustable to various pathologies.



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 42

ART. D700BHE

Claw foot insole, woman template



AVAILABLE COLOURS:


Notes: also available with prefinished base




SIZES: from 34 to 50

ART. 900HEVA

Claw foot insole, man template



AVAILABLE COLOURS:




Notes: also available with prefinished base



SIZES: from 34 to 50

ART. P52001

Insole with a calcaneal shock absorber insert



AVAILABLE COLOURS:




Notes: also available with prefinished base



SIZES: from 34 to 50

ART. P52002/AP

Insole with opening from midfoot to metatarsals for inserting corrective elements



AVAILABLE COLOURS:




Notes:



SIZES: from 34 to 49

ART. COMBI 1

Insole base for the application of corrective elements



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 49

ART. COMBI 2

Insole base for the application of corrective elements



AVAILABLE COLOURS:




Notes: _____



SIZES: from 24 to 46

ART. 38FC12

Composite fiber three-quarter stiff shell



AVAILABLE COLOURS:



Notes: also available with thermoformed metatarsal drop



SIZES: from 34 to 47

ART. PI10/MOUSSE

Pronation wedge: recommended for medial gonalgia



AVAILABLE COLOURS:



Notes: also available with a 6 mm inclined plane




SIZES: from 24 to 46 (only even numbers)

ART. 338

Physiological latex foam insole

Latex foam insole
Density: 300 Kg/m³
Hardness: 25 Shore E

Physiological drop-shaped
metatarsal pad



AVAILABLE COLOURS:



Notes: _____



SIZES: 1st (38/40)
2nd (41/42)
3rd (43/45)


ART. 314

Latex foam insole recommended for metatarsalgia and hammertoes

Latex foam insole
Density: 300 Kg/m³
Hardness: 25 Shore E

Crest pad support

Drop-shaped
metatarsal pad



AVAILABLE COLOURS:



Notes: _____



SIZES: from 34 to 46

ART. COMBI 3

Thin insole recommended for claw foot with metatarsalgia and talalgia treatment



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 46

ART. COMBI 4

Thin insole recommended for claw foot with metatarsalgia and talalgia treatment



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 49

ART. COMBI 5

Thin insole recommended for claw foot with metatarsalgia and talalgia treatment




Notes: _____



SIZES: from 34 to 46

ART. 338X/B

Latex foam comfort insole




Notes: _____



SIZES: from 34 to 47

ART. U8855H

Thin insole recommended for claw foot with metatarsalgia treatment



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 47

ART. U8866H

Thin insole recommended for claw foot with metatarsalgia and talalgia treatment



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 46

ART. 338P

Latex foam comfort insole



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 41

ART. 296A

Woman insole, recommended for mild metatarsalgia



AVAILABLE COLOURS:




Notes: _____




SIZES: from 34 to 46

ART. 338A/B

Latex foam comfort insole



AVAILABLE COLOURS:




Notes: _____



SIZES: from 34 to 46

ART. FLU/SC

Phlebology insole with a calcaneal stabilizer



AVAILABLE COLOURS:



Notes: _____




SIZES: from 34 to 42

COUNTER
DISPLAY



ART. FLD

Phlebology insole



AVAILABLE COLOURS:



Notes: _____




SIZES: from 34 to 47

COUNTER
DISPLAY



ART. LIFE/B

Recommended for backfoot and forefoot pathologies of athlete prevention



AVAILABLE COLOURS:



Notes: _____



SIZES: from 34 to 48

ART. MBDUO

Double density insole recommended for metatarsalgia

Nora[®] Astro Form top layer
Thickness: 6 mm
Density: approx. 200 Kg/m³
Hardness: approx. 8 Shore A

Thermoformed metatarsal bar

Microcel entire base
Thickness: 10 mm
Density: approx. 300 Kg/m³
Hardness: approx. 50 Shore A

AVAILABLE COLOURS:



Notes: _____



SIZES: from 34 to 46

ART. MBE U VARIANT

Recommended for any flat-footedness degrees with mild metatarsalgia treatment

Veolene
Thickness: 23 mm
Density: 140 Kg/m³
Hardness: approx. 29 Shore A

Thermoformed metatarsal bar

AVAILABLE COLOURS:



Notes: also available with backheaded bar or longitudinal arch support only




SIZES: from 34 to 48

ART. MBT BRC VARIANT

Recommended for claw foot with metatarsalgia treatment

Multistratos
Thickness: 15 mm
Density: 180 Kg/m³
Hardness: approx. 35 Shore A

Thermoformed
metatarsal bar



AVAILABLE COLOURS:



Notes: also available with drop-shaped metatarsal pad or longitudinal arch support only



SIZES: from 34 to 50

ART. MBR VLL VARIANT

Recommended for any flat-footedness degrees treatment

Evacolor
Thickness: 22 mm
Density: approx. 200 Kg/m³
Hardness: 35/45 Shore A




AVAILABLE COLOURS:



Notes: also available with drop-shaped metatarsal pad or backheaded bar only

ART. ORTOKIT



Thermoformable material kit for
using foot mold insoles realization

1st size (34 / 37)

2nd size (38 / 41)

3rd size (42 / 46)

Also available with different
combinations of materials

AVAILABLE COLOURS:



SILIPED




Visco-elastic silicone heel cup

Art. 0200
standard

Art. 0250
with a softer central insert

Sizes:
S (33 / 36)
M (37 / 39)
L (40 / 43)
XL (44 / 47)

CALCANEAL WEDGES

**Evacork calcaneal wedges**

Art. R5C Thickness: mm. 5 - Sizes: 1st (35-40) - 2nd (41-46)

Art. R10C Thickness: mm. 10 - Sizes: 1st (35-40) - 2nd (41-46)

Art. R15C Thickness: mm. 15 - Sizes: 1st (35-40) - 2nd (41-46)

Veolform calcaneal wedges

Art. R6C Thickness: mm. 10 - Sizes: 1st (35-40) - 2nd (41-46)

Art. R6CF Thickness: mm. 10 with calcaneal hole

Sizes: 1st (35-40) - 2nd (41-46)

WEDGES TILL METATARSALS

Evacork wedges till metatarsals

Art. R5 Thickness: mm. 5 - Sizes: 1st (24-34) - 2nd (35-40) - 3rd (41-46)

Art. R10 Thickness: mm. 10 - Sizes: 1st (24-34) - 2nd (35-40) - 3rd (41-46)

Art. R15 Thickness: mm. 15 - Sizes: 1st (24-34) - 2nd (35-40) - 3rd (41-46)

Veolform wedges till metatarsals

Art. R65 Thickness: mm. 5 - Sizes: 1st (24-34) - 2nd (35-40) - 3rd (41-46)

Art. R6 Thickness: mm. 10 - Sizes: 1st (24-34) - 2nd (35-40) - 3rd (41-46)

Art. R6F Thickness: mm. 10 with calcaneal hole

Sizes: 1st (24-34) - 2nd (35-40) - 3rd (41-46)



TOTAL INSOLES WITH CALCANEAL WEDGES

**Leather and Evacork total insoles with calcaneal wedges**

Art. RL5/C Thickness mm.5

SIZES: man template (23-47)

woman template (34-42)

Art. RL10/C Thickness mm.10

SIZES: man template (23-47)


woman template (34-42)

Art. RL15/C Thickness mm.15

SIZES: man template (23-47)


woman template (34-42)

ART. 32




Latex foam longitudinal arch support - Sizes: 1st - 2nd

ART. 117




Latex foam longitudinal arch support - Sizes: 1st - 2nd

ART. 189



Latex foam longitudinal arch support - Sizes: 1st - 2nd - 3rd

ART. 30




Latex foam calcaneal support - Sizes: 1st - 2nd

ART. 31



Latex foam calcaneal and longitudinal arch support - Sizes: 1st - 2nd - 3rd

ART. 278




Latex foam calcaneal and longitudinal arch support - Sizes: 1st - 2nd

ART. 185




Latex foam quarter of sphere - Sizes: 1st - 2nd - 3rd

ART. 34



Latex foam crest pad support: Sizes: 1st

ART. 35



Latex foam crest pad support - Sizes: 1st

ART. 113




Latex foam metatarsal pad - Sizes: 1st - 2nd

ART. 114




Latex foam metatarsal pad - Sizes: 1st - 2nd

ART. 427




Latex foam backheaded bar - Sizes: 1st - 2nd - 3rd

ART. 428



Latex foam backheaded bar - Sizes: 1st - 2nd - 3rd

ART. 429



Latex foam backheaded bar - Sizes: 1st - 2nd - 3rd

ART. 2002



Poroten backheaded bar - Sizes: 1st - 2nd - 3rd - 4th

FOOTPRINT DETECTION

**Printfoam**

Footprint detection foam

Size: cm 35x16x5 (box 25 pcs)

Size: cm 31x14x3,5 (box 50 pcs)

Pedifoam

Two-component polyurethane foam (foam + hardener) to realize foot mold

Bottle: Kg. 2

CHECK-UP PAPER

Check up

It allows to realize a fast, accurate and clean podographic examination.

Footprints are printed on an adhesive support and then protected with a transparent sheet. On the back there is a card for the history of pathologies. Box: 50 pcs



GLUE / SOLVENT / ACCESSORIES

**1) Glue**

Glue based on polychloroprene and resins in solution with a medium-long open time. Bottle: Kg. 1,9

2) Glue

Glue based on polychloroprene and resins in solution with a medium-long open time. Bottle: Kg. 16

3) Solvent

Solvent/Thinner for polychloroprenic glue. Bottle: L 1

4) Spalmacol

Glue applicator/spreader with exchangeable brush

5) Glue container

Glue container (L 0,9) with brush

SHAPED COVERING

Shaped covering materials for insoles, available in four different sizes with the possibility to customize them with your logo.

They are realizable with all covering materials in "Materials" section.



1st Size
(23-27)

2nd Size
(28-33)

3rd Size
(34-40)

4th Size
(41-46)

INSOLE BAGS

TNT F/22 easy to use insole bag, available in two colours: light grey and black. We can customize the bag with four-colour printing of your brand.

The maximum printing area is: 10x7 cm - Bag dimensions: 38x16 cm



SA1638



SA1638M


MINIROLL COVERING



- 1) Water-repellent synthetic liner
Colour: Beige
12 cm. x 60 mt
- 2) Water-repellent synthetic liner
Colour: Black
12 cm. x 60 mt
- 3) Alcantara® Podion liner
Colour: Brown
12 cm. x 25 mt
- 4) Nappa synthetic liner
Colour: Ice
12 cm. x 50 mt
- 5) Nappa synthetic liner
Colour: Brown
12 cm. x 50 mt
- 6) Perforated Nappa synthetic liner
Colour: Ice - Black
12 cm. x 50 mt

EVAFORM 100

SINGLE-DENSITY



Colour: Black
 Thickness: 30 mm
 Density: approx. 90 Kg/m³
 Hardness: approx. 25 Shore A
 code 12120 (block cm 35x25)
 code 14120 (block cm 32x25)
 code EVAF100.30N (sheet cm 200x100 approx.)

EVAFORM 120

SINGLE-DENSITY



Colour: Grey
 Thickness: 30 mm
 Density: approx. 120 Kg/m³
 Hardness: approx. 25/30 Shore A
 code 12175 (block cm 35x25)
 code 14175 (block cm 32x25)
 code EVAF120.30 (sheet cm 190x105 approx.)

VEOLENE

SINGLE-DENSITY



Colour: Flesh
 Thickness: 23 mm
 Density: approx. 140 Kg/m³
 Hardness: approx. 29 Shore A
 code 12115 (block cm 35x25)
 code 14115 (block cm 32x25)
 code VN.23C (sheet cm 112x112)

EVAFORM 160


SINGLE-DENSITY



Colour: Beige
 Thickness: 30 mm
 Density: approx. 160 Kg/m³
 Hardness: 30/35 Shore A
 code 12100 (block cm 35x25)
 code 14100 (block cm 32x25)
 code EVAF160.30 (sheet cm 190x105 approx.)

EVAFORM 180


SINGLE-DENSITY



Colour: Light blue
 Thickness: 30 mm
 Density: approx. 180 Kg/m³
 Hardness: 25/35 Shore A
 code 12170 (block cm 35x25)
 code 14170 (block cm 32x25)
 code EVAF180.30 (sheet cm 190x105 approx.)

EVAFORM 240

SINGLE-DENSITY



Colour: Brown
 Thickness: 30 mm
 Density: approx. 240 Kg/m³
 Hardness: 35/40 Shore A
 code 12105 (block cm 35x25)
 code 14105 (block cm 32x25)
 code EVAF240.30 (sheet cm 175x95 approx.)

MICROCEL


SINGLE-DENSITY



Colour: White
Thickness: 25 mm
Density: approx. 300 Kg/m³
Hardness: 50/60 Shore A
code 12110 (block cm 35x25)
code 14110 (block cm 32x25)
code MICRO.25 (sheet cm 160x90 approx.)

EVACOLOR

SINGLE-DENSITY




Colour: Green/Blue marbled
Thickness: 22 mm
Density: approx. 200 Kg/m³
Hardness: approx. 40 Shore A
code 12125 (block cm 35x25)
code 14125 (block cm 32x25)
code EVACOLOR.22 (sheet cm 190x105 approx.)



MULTISTRATOS


SINGLE-DENSITY



Colour: Multicolour
Thickness: 30 mm
Density: approx. 180 Kg/m³
Hardness: app. 35/40 Shore A
code 12140 (block cm 35x25)
code 14140 (block cm 32x25)
code MST.30 (sheet cm 195x100 approx.)

EVAMIX GBC

SINGLE-DENSITY



Colour: Multicolour
Thickness: 30 mm
Density: approx. 180 Kg/m³
Hardness: approx. 25/35 Shore A
code 12140 (block cm 35x25)
code 14140 (block cm 32x25)
code MST.30 (sheet cm 195x100 approx.)

LATEX FOAM

SINGLE-DENSITY



Latex foam
Colour: Grey
Density: approx. 300 Kg/m³
Hardness: approx. 25 Shore E
code LATTICE.15 (block cm 35x25 - thickness 15 mm)
code LATTICE.30 (block cm 35x25 - thickness 30 mm)


EVAFORM 160+240

DOUBLE-DENSITY




BEIGE
Density: approx. 160 Kg/m³
Hardness: 30/35 Shore A
+
BROWN
Density: approx. 240 Kg/m³ - Hardness: 35/40 Shore A
Thickness: 30 mm - Code 13700 (block cm 33x25)

EVAFORM 120+180 DOUBLE-DENSITY



GREY
Density: approx. 120 Kg/m³
Hardness: 25/30 Shore A
+
LIGHT BLUE
Density: approx. 180 Kg/m³ - Hardness: 25/35 Shore A
Thickness: 30 mm - Code 13710 (block cm 33x25)

EVAFORM 180+240 DOUBLE-DENSITY




LIGHT BLUE
Density: approx. 180 Kg/m³
Hardness: 25/35 Shore A
+
BROWN
Density: approx. 240 Kg/m³ - Hardness: 35/40 Shore A
Thickness: 30 mm - Code 13720 (block cm 33x25)

REUMADIAB TRIPLE-DENSITY



FLESH
Thickness: 23 mm
Density: approx. 140 Kg/m³
Hardness: approx. 29 Shore A
+
WHITE
Thickness: 5 mm - Density: approx. 110 Kg/m³ - Hardness: approx. 22 Shore A
+
BLACK - Thickness: 2mm - Density: approx. 180 Kg/m³
Hardness: approx. 40 shore A
code 12300 (block cm 35x25)
code 14300 (block cm 32 x 25)

SPORT TRIPLE-DENSITY




GREEN/BLUE
Thickness: 22 mm
Density: approx. 200 Kg/m³
Hardness: approx. 40 Shore A
+
ANTHRACITE
Thickness: 5 mm - Density: approx. 140 Kg/m³ - Hardness: approx. 29 Shore A
+
LIGHT BLUE - Thickness: 2mm - Density: approx. 180 Kg/m³
Hardness: approx. 40 shore A
code 12305 (block cm 35x25)
code 14305 (block cm 32x25)

HEAVY TRIPLE-DENSITY



WHITE
Thickness: 25 mm
Density: approx. 300 Kg/m³
Hardness: 50/60 Shore A
+
FLESH
Thickness: 5 mm - Density: approx. 140 Kg/m³ - Hardness: approx. 29 Shore A
+
BLACK - Thickness: 2mm - Density: approx. 180 Kg/m³
Hardness: approx. 40 shore A
code 12310 (block cm 35x25)
code 14310 (block cm 32x25)

EVAFORM 160+180+240 TRIPLE-DENSITY




LIGHT BLUE
Density: approx. 180 Kg/m³
Hardness: approx. 25/35 Shore A
+
BEIGE
Density: approx. 160 Kg/m³
Hardness: approx. 30/35 Shore A
+
BROWN - Density: approx. 240 Kg/m³ - Hardness: approx. 35/40 shore A
Thickness: 30 mm
cod. 13730 (block cm. 33x25)

VEOLENE




Polyolefin foam
Density: 140 Kg/m³
Hardness: 29 Shore A
Working temperature: 70° - 90° C
Use: top layer / cover / base
Sheet dimensions: 112 x 112 cm
Thickness: 3/4/15 mm white
1.5/3/4/5/10/23 mm flesh
1.5/4 mm brown
4 mm red
1.5/4 mm green
4 mm blue
3/4/5 mm anthracite

AVAILABLE COLOURS:




PERFORATED VEOLENE



Polyolefin foam
Density: 140 Kg/m³
Hardness: 29 Shore A
Working temperature: 70° - 90° C
Use: top layer / cover
Sheet dimensions: 112 x 112 cm
Thickness: 3/4 mm white
1.5/3/4/5 mm flesh
1.5/4 mm brown
4 mm red
1.5/4 mm green
1.5/4 mm blue
1.5/3/4/5 mm anthracite

AVAILABLE COLOURS:



VEOLFORM



Polyolefin foam
Density: 110 Kg/m³
Hardness: 22 Shore A
Working temperature: 70° - 90° C
Use: Top layer / base
Sheet dimensions: 120 x 120 cm
Thickness: 3 / 5 / 10 mm

PERFORATED VEOLFORM



Polyolefin foam
Density: 110 Kg/m³
Hardness: 22 Shore A
Working temperature: 70° - 90° C
Use: Top layer
Sheet dimensions: 120 x 120 cm
Thickness: 3 / 5 mm

Materials shown on this page and stiff Veolene (page 72) are tested by producers for the manufacture of surface medical devices and, in respect to their biocompatibility, can be expected in compliance with ISO 10993

This has been externally confirmed by independent test institutes on some representative foam grades
-Biological evaluation of medical devices (skin irritation , dermal sensitization) in compliance with ISO 10993
-Resistance against sweat and spittle in compliance with DIN 531600
-Fungi growth resistance in compliance with ISO 846
-SG Approval




EVA 200



EVA Ethylene vinyl acetate
Density: 170/180 Kg/m³
Hardness: 35/45 Shore A
Working temperature: 70° - 90° C
Use: top layer / cover
Sheet dimensions: approx. 190 x 100 cm
Thickness: 1.5mm midnight blue, navy blue, red, black
2mm blue and black
4 mm orange, beige and grey

AVAILABLE COLOURS:




MARBLED EVA 200




EVA Ethylene vinyl acetate
Density: 190/210 Kg/m³
Hardness: 35/45 Shore A
Working temperature: 70° - 90° C
Use: top layer / base
Sheet dimensions: approx. 190 x 100 cm
Thickness: 1.5mm white/light blue
4 mm white/red, black/grey, green/purple,
red/yellow, black/orange
22 mm blue/green and blue/light blue

AVAILABLE COLOURS:




SPOTTED EVA 200



Ethylene vinyl acetate
Density: 190/210 Kg/m³
Hardness: 35/45 Shore A
Working temperature: 70° - 90° C
Use: top layer / cover
Sheet dimensions: approx. 190 x 100 cm
Thickness: 1.5/3/4 mm mimetic
2/4 mm evacolor blue base

AVAILABLE COLOURS:




EVAMIX 180




EVA Ethylene vinyl acetate
Density: 180/200 Kg/m³
Hardness: approx. 35 Shore A
Working temperature: 70° - 90° C
Use: top layer / cover
Sheet dimensions: approx. 190 x 105 cm
Thickness: 1.5/4 mm

GRANITE EVA 160




Ethylene vinyl acetate
Density: approx. 160 Kg/m³
Hardness: approx. 30 Shore A
Working temperature: 70° - 90° C
Use: top layer
Sheet dimensions: approx. 190 x 104 cm
Thickness: 5 mm

NORA® LUNASOFT SL COLOR



EVA Ethylene vinyl acetate
Density: approx. 200 Kg/m³
Hardness: approx. 40 Shore A
Working temperature: 70° - 90° C
Use: top layer
Sheet dimensions: approx. 128 x 89 cm
Thickness: 2 mm orange
4 mm lime

AVAILABLE COLOURS:



NORA® LUNATEC COMBI 4



EVA Ethylene vinyl acetate
Density: approx. 120 (beige) and 230 (white) Kg/m³
Hardness: approx. 22 (beige) and 25 (white) Shore A
Working temperature: 70° - 90° C
Use: top layer / base
Sheet dimensions: approx. 87 x 54 cm
Thickness: 9 mm (3 mm beige + 6 mm white)

NORA® LUNAIRMED




EVA Ethylene vinyl acetate
Density: approx. 80 Kg/m³
Hardness: approx. 16 Shore A
Working temperature: 70° - 90° C
Use: top layer / cover
Sheet dimensions: approx. 92 x 61 cm
Thickness: 3/4 mm

NORA® LUNALASTIK



EVA Ethylene vinyl acetate
Density: approx. 230 Kg/m³
Hardness: approx. 25 Shore A
Working temperature: 70° - 90° C
Use: top layer
Sheet dimensions: approx. 85 x 61 cm
Thickness: 4 mm

NORA® LUNASOFT SL TRENDLINE



EVA Ethylene vinyl acetate
Density: approx. 200 Kg/m³
Hardness: approx. 40 Shore A
Working temperature: 70° - 90° C
Use: cover
Sheet dimensions: approx. 128 x 89 cm
Thickness: 2 mm

NORA® LUNATEC EP



EVA Ethylene vinyl acetate
Density: approx. 200 Kg/m³
Hardness: approx. 22 Shore A
Working temperature: 70° - 90° C
Use: cover
Sheet dimensions: approx. 88 x 59 cm
Thickness: 2 mm


NORA® LUNATUR 18 WALNUT



EVA Ethylene vinyl acetate + walnut shells (more than 20%)
Density: approx. 120 Kg/m³
Hardness: approx. 18 Shore A
Working temperature: 70° - 90° C
Use: top layer / cover
Sheet dimensions: approx. 98 x 64 cm
Thickness: 3 mm



NORA® LUNASOFT SL COLOR PLUS




EVA Ethylene vinyl acetate
Density: approx. 200 Kg/m³
Hardness: approx. 40 Shore A
Working temperature: 70° - 90° C
Use: cover
Sheet dimensions: approx. 128 x 89 cm
Thickness: 2 mm

EVAFLEX



EVA Ethylene vinyl acetate
Density: approx. 230 Kg/m³
Hardness: approx. 35 Shore A
Working temperature: 70° - 90° C
Use: cover
Sheet dimensions: approx. 175 x 95 cm
Thickness: 1.5/2 mm marbled beige/brown
1.5 mm marbled yellow/grey


AVAILABLE COLOURS:



BM

GG

NORA® ASTRO FORM 8



Light cellular rubber
Density: approx. 200 Kg/m³
Working temperature: 70° - 90° C
Use: top layer / padding
Sheet dimensions: approx. 77 x 56 cm
Thickness: 6 mm

NORA® ASTRO MED 10




Closed cell rubber
Density: approx. 270 Kg/m³
Hardness: approx. 10 Shore A
Working temperature: not thermoformable
Use: cover / padding
Sheet dimensions: approx. 105 x 70 cm
Thickness: 3 mm

PORON®




Urethane foam
Density: approx. 270 Kg/m³
Hardness: approx. 20 Shore 0
Working temperature: not thermoformable
Use: cover / padding
Sheet dimensions: approx. 137 x 73 cm
Thickness: 1.6 mm

RED MEMORY



Urethane foam
Density: approx. 240 Kg/m³
Hardness: approx. 18 Shore 0
Working temperature: not thermoformable
Use: cover / padding
Sheet dimensions: approx. 137 x 73 cm
Thickness: 3 / 6 mm

TRANSPIRA + MEMORY



Modified SBR foam
Density: approx. 230 (light blue) + 180 (red) Kg/m³
Working temperature: not thermoformable
Use: top layer / padding
Sheet dimensions: approx. 100 x 100 cm
Thickness: 4.5 mm

AIRLASTIC



Urethane foam
Density: approx. 320 Kg/m³
Hardness: approx. 17 Shore 0
Working temperature: not thermoformable
Use: cover / padding
Sheet dimensions: approx. 137 x 73 cm
Thickness: 1.5/3/6 mm

PERFORATED AIRLASTIC



Urethane foam
Density: approx. 320 Kg/m³
Hardness: approx. 17 Shore 0
Working temperature: not thermoformable
Use: cover / padding
Sheet dimensions: approx. 137 x 73 cm
Thickness: 1.5/3 mm

COVERED PORON®



Urethane foam + Microfibre
Density: approx. 270 Kg/m³
Hardness: approx. 20 Shore 0
Working temperature: not thermoformable
Use: cover
Sheet dimensions: approx. 137 x 73 cm
Thickness: 2,1 mm

COVERED AIRLASTIC




Urethane foam + Microfibre
Density: approx. 320 Kg/m³
Hardness: approx. 17 Shore 0
Working temperature: not thermoformable
Use: cover
Sheet dimensions: approx. 137 x 73 cm
Thickness: 3,7 mm

PLASTAZOTE



Polyethylene foam
Density: approx. 45 Kg/m³
Hardness: approx. 62 Shore 00
Working temperature: 60° - 70° C
Use: top layer / cover
Sheet dimensions: approx. 200x100 cm
Thickness: 3/5/10 mm white
3 mm self-adhesive white
3 mm black
4 mm red
4 mm yellow

AVAILABLE COLOURS:



PERFORATED PLASTAZOTE



Polyethylene foam
Density: approx. 45 Kg/m³
Hardness: approx. 62 Shore 00
Working temperature: 60° - 70° C
Use: top layer / cover
Sheet dimensions: approx. 200 x 100 cm
Thickness: 3/5 mm

STIFF VEOLENE



Polyolefin foam
Density: 170 Kg/m³
Hardness: 40 Shore A
Working temperature: 80° - 100° C
Use: base
Sheet dimensions: 100x100 cm
Thickness: 5/7/10/20 mm

MICROCEL




EVA Ethylene vinyl acetate
Density: approx. 300 Kg/m³
Hardness: approx. 50 Shore A
Working temperature: 80° - 100° C
Use: base
Sheet dimensions: approx. 160 x 90 cm
Thickness: 5/10/15/25 mm white
5/10 mm orange
10 mm light grey

AVAILABLE COLOURS:



DUROCEL




EVA Ethylene vinyl acetate
Density: approx. 320 Kg/m³
Hardness: approx. 60 Shore A
Working temperature: 80° - 100° C
Use: base
Sheet dimensions: approx. 165 x 90 cm
Thickness: 3/5/10 mm black
3 mm yellow, red
5 mm light green

AVAILABLE COLOURS:




MULTISTRATOS



EVA Ethylene vinyl acetate
Density: approx. 180 Kg/m³
Hardness: approx. 35 Shore A
Working temperature: 80° - 100° C
Use: base
Sheet dimensions: approx. 195 x 105 cm
Thickness: 15/30 mm

MULTISTRATOS NATURE




EVA Ethylene vinyl acetate
Density: approx. 180 Kg/m³
Hardness: approx. 35 Shore A
Working temperature: 80° - 100° C
Use: base
Sheet dimensions: approx. 195 x 105 cm
Thickness: 12 mm

EVAFORM 160



EVA Ethylene vinyl acetate
Density: approx. 160 Kg/m³
Hardness: approx. 30 Shore A
Working temperature: 70° - 90° C
Use: base
Sheet dimensions: approx. 190 x 105 cm
Thickness: 10/30 mm

EVAFORM 240



EVA Ethylene vinyl acetate
Density: approx. 240 Kg/m³
Hardness: approx. 35/40 Shore A
Working temperature: 70° - 90° C
Use: base
Sheet dimensions: approx. 190 x 105 cm
Thickness: 10/30 mm

EVACORK




EVA Ethylene vinyl acetate
Density: approx. 240 Kg/m³
Hardness: approx. 50 Shore A
Working temperature: 70° - 90° C
Use: base
Sheet dimensions: approx. 175 x 115 cm
Thickness: 3/5/7/10/15 mm

SUPERSUGHERO



EVA Ethylene vinyl acetate
Density: 260 Kg/m³
Hardness: 60 Shore A
Working temperature: 70° - 90° C
Use: base
Sheet dimensions: approx. 175 x 95 cm
Thickness: 5/10 mm pink
7/10 mm light blue

AVAILABLE COLOURS:



SUGHERO 839D



EVA Ethylene vinyl acetate
Density: approx. 440 Kg/m³
Hardness: approx. 72 Shore A
Working temperature: 70° - 90° C
Use: base
Sheet dimensions: approx. 90 x 70 cm
Thickness: 5/10 mm

CORK MOUSSE



Natural rubber foam and cork
Density: approx. 420 Kg/m³
Hardness: 20/35 Shore A
Working temperature: not thermoformable
Use: top layer / base
Sheet dimensions: 90 x 60 cm
Thickness: 3/6/10 mm

THERMO MOUSSE



Natural rubber and thermoplastic copolymers mixture
Density: approx. 490 Kg/m³
Working temperature: 90° - 100° C
Use: top layer / base
Sheet dimensions: 90x60cm
Thickness: 3/6/10 mm

POROTEN



Foam rubber
Density: approx. 130 Kg/m³
Hardness: approx. 35 Shore 00
Working temperature: not thermoformable
Use: padding
Sheet dimensions: approx. 200x100cm
Thickness: 5 / 10 mm

PEDIFLEX



Resin on polyester fiber fabric
Working temperature: 70° - 80° C
Use: top layer / base
Sheet dimensions: 150 x 100 cm
Thickness: 0,8 mm standard
0,8 mm perforated

TALIN



Resin on polyester fiber fabric
Weight: approx. 730 g/m²
Working temperature: 70° - 80° C
Use: top layer / base
Sheet dimensions: 150x100 cm
Thickness: approx. 0,85 mm

STARFLEX



Resin on polyester fiber fabric
Weight: approx. 510 g/m²
Working temperature: 70° - 80° C
Use: bottom layer covering
Sheet dimensions: 150x100 cm
Thickness: approx. 0,65 mm

SINTEX (one glued side)



Resin on polyester fiber fabric
Working temperature: 80° - 90° C
Use: reinforcement
Sheet dimensions: 150 x 100 cm
Thickness: approx. 1,0 / 1,4 / 1,8 mm

SINTEX (two glued sides)



Resin on polyester fiber fabric
Working temperature: 80° - 90° C
Use: reinforcement
Sheet dimensions: 150 x 100 cm
Thickness: approx. 1,0 / 1,4 mm

DUROFORM




Resin on polyester fiber fabric
Working temperature: 80° - 90° C
Use: reinforcement
Sheet dimensions: 150 x 100 cm
Thickness: approx. 0,8 / 1,2 / 1,4 mm

RESINPLUS



Counter-material on thermoplastic polyester base
Working temperature: 80° - 90° C
Use: reinforcement
Sheet dimensions: 150 x 100 cm
Thickness: approx. 2,1 mm white
approx. 1,2 mm magenta
approx. 1,0 mm green

AVAILABLE COLOURS:



NOENE®



Synthetic cellular rubber with TNT fabric
Density: approx. 640 Kg/m³
Hardness: approx. 34 Shore A
Working temperature: not thermoformable
Use: Shock Absorber
Sheet dimensions: approx. 130x77 cm
Thickness: 2 mm

NORA® AERO SORB M



Foam cellular rubber
Density: approx. 160 Kg/m³
Working temperature: 70° - 90° C
Use: Shock absorber
Sheet dimensions: approx. 80 x 55 cm
Thickness: 2 mm

PORON® XRD



Urethane foam
Density: approx. 240 Kg/m³
Hardness: approx. 32 Shore 0
Working temperature: not thermoformable
Use: Shock absorber
Sheet dimensions: approx. 137 x 73 cm
Thickness: 3 mm

PORON® ETH VIVE



Open-cell urethane foam
Density: approx. 256 Kg/m³
Working temperature: not thermoformable
Use: Cover / Padding
Sheet dimensions: approx. 137x730 cm
Thickness: 3 mm

POR 251



Closed cell rubber
Density: approx. 250 Kg/m³
Hardness: 65/70 Shore 00
Working temperature: 70° - 90° C
Use: Shock absorber
Sheet dimensions: approx. 200x100 cm
Thickness: 2 mm

LATEX FOAM



Synthetic and natural latex mixture
Density: approx. 300 Kg/m³
Hardness: approx. 25 Shore E
Working temperature: not thermoformable
Use: top layer / base
Sheet dimensions: approx. 44x22 cm
Thickness: 6 / 10 / 15 mm

LATEX



Latex on cotton fabric
Working temperature: not thermoformable
Use: soft padding
Sheet dimensions: approx. 150 x 100 cm
Thickness: 6 / 10 mm

GOATSKIN



Vegetable tanned goatskin
Use: cover
Leather dimension: approx. 0,5 m²
Thickness: approx. 0,5 mm

PERFORATED GOATSKIN



Perforated vegetable tanned goatskin
Use: cover
Leather dimension: approx. 0,5 m²
Thickness: approx. 0,5 mm

CALFSKIN



Vegetable tanned calfskin
Use: cover
Leather dimension: approx. 1,2 m2
Thickness: approx. 0,9 mm

VEGETABLE TANNED SUEDE



Vegetable tanned suede
Use: cover
Leather dimension: approx. 0,5 m2
Thickness: approx. 0,5 mm

PIGSKIN



Vegetable tanned pigskin
Use: Cover
Leather dimension: approx. 1,4 m2
Thickness: approx. 0,8mm

LEATHER




First choice vegetable tanned leather
Use: Top Layer
Leather dimension: approx. 1,2 m2
Thickness: approx. 1,1 / 1,5 / 1,8 / 2,5 mm

ALCANTARA® PODION



AVAILABLE COLOURS:




Polyester and polyurethane microfiber
Weight: approx. 230 g/m2
Use: Cover
Sheet dimensions: approx. 145x100 cm
Thickness: approx. 0,8mm

SUEDE



AVAILABLE COLOURS:




Nylon and polyurethane microfiber
Weight: approx. 270 g/m2
Use: Cover
Sheet dimensions: approx. 145x100 cm
Thickness: approx. 0,8 mm

DELIKATESSE VACCHETTA LINING



Synthetic lining made of viscose and polyurethane
Weight: approx. 290 g/m²
Use: Cover
Sheet dimensions: approx. 145x100 cm
Thickness: approx. 0,75 mm

CARBON LINING




Synthetic lining made of cotton and polyurethane
Weight: approx. 325 g/m²
Sheet dimensions: approx. 140x100 cm
Thickness: approx. 0,8 mm

NAPPA LINING



Synthetic lining made of viscose and polyurethane
Weight: approx. 225 g/m²
Use: Cover
Sheet dimensions: approx. 145x100 cm
Thickness: approx. 0,75 mm

AVAILABLE COLOURS:




WATER-REPELLENT LINING



Synthetic lining made of cotton, polyester and polyurethane
Weight: approx. 290 g/m²
Sheet dimensions: approx. 140x100 cm
Thickness: approx. 0,8 mm

AVAILABLE COLOURS:



PRESTIGE LINING




Synthetic lining made of viscose and polyurethane
Weight: approx. 240 g/m²
Use: Cover
Sheet dimensions: app. 145x100 cm
Thickness: approx. 0,75 mm

DELUXE LINING



Synthetic lining made of viscose and polyurethane
Weight: approx. 230 g/m²
Use: Cover
Sheet dimensions: approx. 145x100 cm
Thickness: approx. 0,75 mm

COMPOSITE FIBER



Thermoplastic polyurethane and glass fiber
Weight: approx. 2250 g/m²
Working temperature: 200° C for max 4 minutes
Use: top layer / reinforcement
Sheet dimensions: 100x62 cm
Thickness: approx. 1,25 mm

POLYETHYLENE



Extruded polyethylene
Density: approx. 950 Kg/m³
Hardness: approx. 62 Shore D
Working temperature: 170° / 190° C
Use: Wraparound shells / orthopedic corset
Sheet dimensions: 200x100 cm
Thickness: 2/3/4/5 mm

SPONGE RUBBER




Natural rubber
Density: approx. 700 Kg/m³
Working temperature: not thermoformable
Use: Corrections/supports
Sheet dimensions: 100x50 cm
Thickness: 3/5 mm

SILICONE RUBBER



Platinic silicone sheet high breaking strenght
Weight: approx. 2250 g/m²
Use: Silicone rubber for vacuum
Sheet dimensions: on customer request
Thickness: 2 mm

(Mod. V7-2 rev. 0)



Date:

Address:

Delivery address:

Tel

Fax

VAT ID No.

Bank:

Payment:

Mail

INSOLES

[illegible][illegible][illegible]

SHAPED COVERING	23-27	28-33	34-40	41-46

CUSTOMER SIGNATURE / STAMP

GENERAL TERMS OF SALE: In case of failed or incomplete indication of item code, the item with closest description will be sent. Order will be automatically confirmed by Pedsan if the customer doesn't receive other information within 5 days of receipt of the order. Ordered items will be delivered within 30 days from the order date, except in cases of force majeure. The lack of items or the delay in delivery not entitle the customer to apply for refunds or payment suspensions. Goods are shipped with free specific packaging. Any claims must be made within 8 days of receipt of the goods. Any returns must be previously authorized by our head office and always shipped back in free port with delivery note. The Court of Perugia shall have exclusive jurisdiction over any disputes relating to the terms and conditions mentioned above.

Lined area for notes.

QUALITY

All insoles and dysmetria wedges manufactured by Pedsan s.r.l. are certified CE Class 1 according to the Legislative Decree no. 46 of 24th February 1997 (implementation of Directive 93/42 EEC)

GENERAL TERMS OF SALE: In case of failed or incomplete indication of item code, the item with closest description will be sent. Order will be automatically confirmed by Pedsan if the customer doesn't receive other information within 5 days of receipt of the order. Ordered items will be delivered within 30 days from the order date, except in cases of force majeure. The lack of items or the delay in delivery not entitle the customer to apply for refunds or payment suspensions. Goods are shipped with free specific packaging. Any claims must be made within 8 days of receipt of the goods. Any returns must be previously authorized by our head office and always shipped back in free port with delivery note. The Court of Perugia shall have exclusive jurisdiction over any disputes relating to the terms and conditions mentioned above.

NOTICE

The information in this catalogue have the sole purpose of illustrating products and they are intended as mere description of products which do not guarantee their features.

The images published in this catalogue are purely indicative: colours may be different due to the printing process.

©COPYRIGHT PEDSAN

All rights reserved. Any partial or entire use of contents of this catalogue is prohibited, including storage, reproduction, revision, dissemination or distribution of contents through any technology platform, support or network, without prior authorisation by Pedsan srl.

Pedsan

PEDSAN s.r.l.
Via Sandro Penna, 112
06132 - S. Andrea delle Fratte
PERUGIA - Italy
Tel. +39 075.5289118
Fax +39 075.5271853
www.pedsan.it
www.pedsan.com
info@pedsan.it

